

Case Share

- Click the patient's name from the dashboard or find them via search patients
- Click the "+" in the Case Sharing section

Case Sharing



Case has not been shared.

Case Share is a HIPAA compliant way to share your patient's BeamReaders profile. When you use Case Share you'll give access to all the patient files (DICOM, 2D images, etc) and any BeamReaders services (radiology report, image portfolio). You can also use Case Comments to have a secure conversation about the patient and their treatment options. Case Share is meant only for doctor-to-doctor sharing, not patient access.

Don't worry if you want to share with someone that doesn't already have a BeamReaders account they will receive a personalized invite and will be able to setup a free account.

Case Share doesn't cost anything if included with a BeamReaders service otherwise it is only a buck!

STEP 1

Choose Connections




Share this Case

Test Customer
#300011

- 1 Choose Connections**
Your practice connections
- 2 Add a Connection**
Add to your connections
- 3 Add a Note**
- 4 Review**

The next step allows you to add a new connection if the person you would like to share with isn't in the list below.

Search connections

Practice Connections		
	Bob Parr affirmaBeamer2@gmail.com Practice 1	<input checked="" type="checkbox"/>
	Dash Parr affirmaBeamer4@gmail.com Practice 1	<input type="checkbox"/>
	Violet Parr affirmaBeamer2@gmail.com	<input type="checkbox"/>

CANCEL NEXT

→ You'll build a list of connections from prior shares. If the person you want to share with is already in your connections list select them

STEP 2

Add a Connection

Share this Case
Test Customer #300011

✓ Choose Connections
Your practice connections

2 Add a Connection
Add to your connections

3 Add a Note

4 Review

Optional: create a new connection to share this case with.

<small>Title</small>	<small>First Name *</small>	<small>Last Name *</small>
	▼ New	Connection
<small>Email</small>		
newconnection@beamreaders.com		

*If you need to add additional new connections to this case you can repeat this work-flow.

CANCEL
PREVIOUS
NEXT

→ To add a new doctor to share with all you need is to enter their name and e-mail

STEP 3

Add a Note

Share this Case

Test Customer
#3000011

✓ Choose Connections — ✓ Add a Connection — 3 Add a Note — 4 Review

Optional: Add a personal note to the share invitation e-mail. This note will not be part of the case record.

Thanks for the referral! Here is the patient's 3D scan and a great report from BeamReaders that unlocks the mysteries of the scan.

CANCEL PREVIOUS NEXT

→ Add an introductory comment to give context to those you're sharing with

STEP 4

Review

Share this Case

Test Customer
#300011

✓ Choose Connections — ✓ Add a Connection — ✓ Add a Note — 4 Review

Your case will be shared with the following people:

Bob Parr (affirmaBeamer2@gmail.com)
New Connection (newconnection@beamreaders.com)

Thanks for the referral! Here is the patient's 3D scan and a great report from BeamReaders that unlocks the mysteries of the scan.

Click "Share" to send your sharing invitations. These users will have access to view your case, respond to comments, and download any files.

CANCEL PREVIOUS SHARE

- Double check who you're sharing with and your note
- Click Share!
 - Everyone will receive an e-mail notifying them about the share.
 - If the doctor already has a BeamReaders account the case will show up on their dashboard under Shared Cases.
 - If the doctor is new to BeamReaders they will receive an email introducing them to BeamReaders and helping them setup their account so they can securely access the patient's profile.